

# maxill H<sub>2</sub>O Troubleshooting Guide

Observation	Result	Action code	Steps
< 500 CFU/mL	PASS	1	<ul style="list-style-type: none"> <li>• Continue with best practices</li> </ul>
500 - 1000 CFU/mL	CAUTION	2	<ul style="list-style-type: none"> <li>• Perform a pre-clean with Zymax</li> <li>• Perform a disinfect with Aqua-FX</li> <li>• Clean and disinfect your DUWL and water bottle</li> </ul>
> 1000 CFU/mL	FAIL	3	<ul style="list-style-type: none"> <li>• Perform a pre-clean with Zymax and a disinfect with Aqua-FX. Multiple treatments may be necessary</li> <li>• Consider source water testing and dental water line unit maintenance activities</li> <li>• Ensure water sampling technique is not contaminating samples</li> <li>• Confirm effective cleaning / PAA treatments through additional testing</li> </ul> <p>CAUTION: If repeated observations of high bacterial counts occur when using best practices, consult with your DUWL service technician</p>

## maxill DUWL Treatment Kit

ITEM #81325

The maxill DUWL Treatment Kit was designed with both cleaning and best practice maintenance in mind. This weekly, 2 Step system includes a Zymax pre-cleaning step to erode biofilms and an Aqua-FX disinfection to oxidize bacteria in the DUWLs. Zymax should be used daily for best practice maintenance. Use in combination with maxill DUWL in-office or lab testing options for a complete end-to-end monitoring and maintenance solution.

**NO** Tablets      **NO** Silver      **NO** Reverse Osmosis Units  
**NO** Chemicals      **NO** Zinc      **NO** UV Lighting Units  
**NO** Filters      **NO** Iodine      **NO** Inline Disks

### Kit Includes:

- Zymax 128 Enzymatic Cleaning Solution
- Zymax Chairside Bottle
- maxill Aqua-FX
- maxill Aqua-FX Measuring Cup
- Peracetic Acid Test Strips
- Aqua-FX Chairside Bottle

**maxill H<sub>2</sub>O = WATER & ONLY WATER**

