

ECONOMICAL 64x CONCENTRATE

Evaculine (6) **DENTAL EVACUATION SYSTEM CLEANER**

- Fresh lemon scent
- Non-acid
- Removes organic waste
 Cleans and deodorizes
 - Corrosion protection
 - Low foaming formulation

When diluted, *Evaculine 64* is designed to clean and deodorize dental institutions' evacuation systems. Evaculine 64 is a non-acid concentrate to be diluted with water.

Directions: Dilute 1 part *Evaculine 64* with 63 parts water, near room temperature. E.g. Mix 15 mL $/\frac{1}{2}$ fl. oz. of Evaculine 64 with water making up a working solution to 0.95 L / 1 gt. Aspirate or vacuum the mixed solution through the evacuator lines and cuspidor. Repeat this procedure for each dental unit at least once daily.

1 L 33.8 FL. OZ. (1.06 QT.)

ITEM 60036

Evaculine 64 - Dental Evacuation System Cleaner

DANGER

Causes severe skin burns and eye damage. May cause an allergic skin reaction.





PREVENTION

Do not breathe mist or vapour.

Contaminated work clothing should not be allowed out of the workplace.

Wash thoroughly after handling.

Wear protective gloves/protective clothing/eye protection/face protection.

RESPONSE

IF ON SKIN (or hair): Take off immediately all contaminated clothing.

Rinse skin with water.

Immediately call a POISON CENTER/doctor.

If skin irritation or rash occurs: Get medical advice/attention.

Wash contaminated dothing before reuse.

Specific treatment (see information on this label).

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

IF SWALLOWED: Rinse mouth, Do NOT induce vomiting.

STORAGE

Store locked up.

DISPOSAL

Dispose of contents/container in accordance with local/regional/national/international regulations.

Made in Canada with domestic and imported parts



Manufactured for: | USA: maxill inc. Cortland, OH 44410 USA

Canada: maxill inc. St. Thomas, ON N5R 6C8 Canada





+D89260705600360K