Onyx™

INSTRUCTIONS FOR USE

40% Phosphoric Acid Etching Gel

GB - ENGLISH

INSTRUCTIONS FOR USE

Product Description: Onyx™ is a viscous, high-contrast 40% phosphoric acid etchant (H3P04). The black silicon carbide particles provide easier visiblity for placing the etchant and then helps verify that the etchant has been cleanly and completely rinsed away.

1. Prepare tooth using pumice, disc or diamond bur on uncut enamel. Rinse and dry thoroughly.
2. Remove cap from syringe and screw on a disposable needle tip.
3. Apply Onyx etching gel to enamel and/or dentin.
4. Let stand on tooth surface for 15 seconds. Rinse with water and high-volume suction.
5. Dry gently. Tooth surface should appear dull. If not, repeat above steps.
6. Remove needle tip, re-cap and disinfect syringe barrel.

INDICATIONS: Appropriate for any bonding system requiring selective or total etching.

CAUTION & WARNINGS:

- Onyx contains 40% phosphoric acid; if in contact with any soft tissue, IMMEDIATELY rinse area with copious amounts of water.
- Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check.
- Re-cap, disinfect and wipe syringe with an intermediate level disinfectant between uses.
- Dispose of tip and empty syringe properly.
- Near pulp exposures should be treated with a protective base prior to placing etchant.
- Do not use after expiration date noted on container.

DISCLAIMER: Centrix products have been developed for use in dentistry. We have provided written, verbal and online information as to the use of each product. This information does not relieve you from examining our products upon receipt to determine their suitability for use in the intended manner of application. The usage of our products is beyond our control and the user maintains full responsibility. We guarantee the quality of our products in accordance with the existing standards and corresponding to the conditions stated in our instructions for use and safety data sheets.

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SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Trade Name: Onyx, Etching Gel
Synonyms: 40% Phosphoric Acid etching gel
Part (Item) Number: See Above
Application: Etching gel for dental restorations
Manufacturer: Centrix Inc.
Address: 770 River Road
City, State, Zip: Shelton Ct, 06484
Emergency Telephone Number: 203-929-5582
Information Telephone Number: 800-235-5862

SECTION II – HAZARD IDENTIFICATION
Classification of the Substance or Mixture:
Acute toxicity, Oral (Category 4)
Acute toxicity, Dermal (Category 5)
Serious eye damage (Category 1)

Risk Phrases:
Symbol: C
R34: Causes burns.
Label Elements:
Trade Name: PHOSPHORIC ACID
Signal Word: Danger

Hazard Statements:
H302: Harmful if swallowed.
H313: May be harmful in contact with skin.
H314: Causes severe skin burns and eye damage.
H330: Fatal if inhaled.
Precautionary Statements:
P260: Do not breathe dust / fume / gas / mist / vapors / spray.
P280: Wear protective gloves / protective clothing / eye protection / face protection.
P284: Wear respiratory protection.
P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing.
P310: Immediately call a POISON CENTER or doctor / physician

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS
Hazardous Components: Phosphoric Acid
Concentration: 40%
CAS No.: 7664-38-2
EC Number: 231-633-2
Index Number: 015-011-00-6
Molecular Weight: 98 g/mol

SECTION IV - FIRST AID MEASURES

Skin: Immediately wash skin with plenty of water.
If irritation occurs seek medical advice.
Ingestion: If swallowed do not induce vomiting. Give large quantities of water and milk.
Seek medical attention.
Eyes: Flush with plenty of water.
Seek immediate medical attention by an Ophthalmologist.

SECTION V - FIRE-FIGHTING MEASURES

The product is not considered to be a fire hazard. The information below are given for extinguishing a fire generated by other sources.
Suitable Extinguishing Media: Water, Chemical Foam, Carbon Dioxide, or Dry Chemical
Forbidden Extinguishing Media: None.
Special Exposure Hazards: Toxic fumes of Phosphoric Acid.
Unusual fire and explosion hazards: Contact with some metals may cause formation of hydrogen gas.
Special Protection Equipment: Self Contained Breathing Apparatus if toxic fumes are generated.

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal Precautions: See sections VII and VIII.
Environmental Precautions: None.
Methods for Cleaning up: Clean up using impervious gloves and dispose of in an approved manner.

SECTION VII - HANDLING AND STORAGE

Handling: Avoid unnecessary exposure.
Storage: Room Temperature.
Containers: Supplied by manufacturer.
Incompatibility: Strong caustics, sulfites, sulfides and metals.
Other Precautions: Use according to directions.
Avoid contact with eyes, skin and soft tissue.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION
General Measures Ventilation: N/A
Limits: TWA: 1mg/m3; TLV: 3 mg/m3.
Respiratory Protection: N/A
Eye Protection: Safety Glasses or Goggles.
Skin Protection: Respect good personal hygiene practice.
Avoid unnecessary exposure.
Hands Protection: Impervious gloves
Other Equipment: Lab Coat.

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Gel
Explosion Limits: N/A
Odor: None
pH 10% solution: below 2
Boiling Point: 100°C
Flash Point: N/A

SECTION X - STABILITY AND REACTIVITY.

Stability: Stable.
Conditions to avoid: Temperatures over 40°C
Materials to avoid: Strong Caustics, Sulfites, Sulfides, Metals.
Hazardous Decomposition Products: Toxic fumes if temperature is above 100°C.
Hazardous Polymerization Products: None

SECTION XI - TOXICOLOGICAL INFORMATION.

Ingestion: May cause sore throat, abdominal pain, nausea and severe burns of mouth, throat and stomach.
Eye contact: May cause redness, pain, blurred vision, eye burns or permanent eye damage.
Skin Contact: May cause redness, pain and severe skin burns.
Inhalation: Inhalation is not a hazard unless the product is heated at high temperatures.
In these conditions may cause irritation of the respiratory system, asphyxia or permanent damage to the lungs.
Prolonged Exposure: Burns to eye, soft tissue or respiratory tract.

SECTION XII - ECOLOGICAL INFORMATION.

There is no specific information on the product. Use according to good working practice. Avoid releasing the product into environment.

SECTION XIII - DISPOSAL CONSIDERATIONS.
Dispose in accordance with local authorities.

SECTION XIV - TRANSPORT INFORMATION.
Substances determining Hazard: Phosphoric Acid.
SEA TRANSPORTATION: IMCO # 8204 UN # 1805
AIR TRANSPORTATION: ICAO/IATA CLASS: 8
GROUND TRANSPORTATION: RID/ADR CLASS: 9
SUBSTANCE IDENTIFICATION: # 1805 ITEM # 11C
HAZARD LABEL: 9
HAZARD IDENTIFICATION: 80

SECTION XV - REGULATORY INFORMATION
Hazard Classification Symbol: C
R34 Causes burns
S24/S25 Avoid contact with skin and eyes
S26/S28 In case of contact with eyes or skin, rinse immediately with plenty of water and seek medical advice.

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION:

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

centrix®
MAKING DENTISTRY EASIER

	CENTRIX INCORPORATED 770 River Road, Shelton, CT 06484 USA (203) 929-5582 • (800) 235-5862 • FAX (203) 929-6804 www.centrixdental.com
	
	Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden
	
	EMERGO AUSTRALIA Level 20, Tower II, Darling Park, 201 Sussex Street Sydney, NSW 2000, Australia