maxill

DIN 02488973

Minuteman NEXGEN

DISINFECTANT WIPES

MYCOBACTERICIDAL • VIRUCIDAL¹ • BACTERICIDAL • FUNGICIDAL

Disinfectant-Sanitizer Intermediate Level Hospital Disinfectant Wipes

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n-Alkyl (50% C14, 40% C12, 10% C16)	
dimethyl benzyl ammonium chloride...	0.34% w/w
Didecyl dimethyl ammonium chloride	0.51% w/w
Ethanol	6.84% w/w
Isopropanol	13.68% w/w

160

WIPES/
LINGETTES

15 cm x 18 cm / 6" x 7"

WARNING

PRECAUTIONARY STATEMENTS:

Hazards to human and domestic animals: Flammable liquid and vapour. Avoid breathing mist or vapour. Harmful in contact with skin. Causes skin irritation. Causes serious eye irritation. May cause drowsiness or dizziness. Use only outdoors or in a well-ventilated area. Wear protective gloves and eye protection. Keep away from heat and open flames. No smoking. Wash hands thoroughly after handling.

FIRST AID:

IF ON SKIN (or hair): Wash with plenty of water. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention. Call a POISON CENTER/doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Get medical help if you feel unwell. Take the container label or product name and Drug Identification Number (DIN) with you when seeking medical attention.

STORAGE: Store in a well-ventilated place. Keep cool. Store locked up. Keep container tightly closed. Do not reuse or refill container except as described in the directions for use. (Refill only with this product.)

DISPOSAL: Recycle this container or discard container in trash. Do not reuse wipe. Dispose of wipe in trash after use. Do not flush in the toilet. For information on disposal of unused, unwanted product, contact the manufacturer or the appropriate municipal/provincial/territorial agency or product stewardship organization.

¹The empty canister from product #61160 is reusable when refilled with product labeled as refill pouch, #61190.

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#61160

KEEP OUT OF REACH OF CHILDREN
READ THE LABEL BEFORE USING

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Made in Canada with domestic and imported ingredients

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[REF] 6070561160

D89260705611600J

tb Minuteman NEX GEN Disinfectant Wipes is intended for use on pre-cleaned, hard non-porous inanimate surfaces in areas where the disinfection of surfaces is required. When used as directed tb Minuteman NEX GEN Disinfectant Wipes cleans and disinfects hard, non-porous, non-food contact surfaces found in ambulances, dental facilities, hospitals, nursing homes, paramedic units, schools, clinics, parlors, bathrooms or any other places that are concerned with the hazards of cross-contamination from infectious microorganisms. When used as directed tb Minuteman NEX GEN Disinfectant Wipes may be used on the following hard, non-porous, non-food contact surfaces: bed & security rails, cabinets, commodes, counters, equipment surfaces & knobs, pet care surfaces, tables, toys, vehicle interior surfaces. When used as directed tb Minuteman NEX GEN Disinfectant Wipes may be used on the following hard, non-porous, non-critical parts of the following items/surfaces: ambulance equipment, blood pressure devices, CPR training devices, medical equipment and dental chairs¹, EMS equipment, exam tables, respiratory care & therapy equipment and any surfaces that are potentially contaminated with disease-causing microorganisms. If you are uncertain as to the material composition of your surface, confirm compatibility with the manufacturer before proceeding. For plastic and painted surfaces, spot test on an inconspicuous area before use.

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1ORGANISMS CONTROLLED / CONTACT TIMES: When used as directed on pre-cleaned hard, non-porous non-food contact surfaces, tb Minuteman NEX GEN Disinfectant Wipes is a sanitizer and disinfectant, effective against the following microorganisms.

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Rotavirus
Human Coronavirus
SARS-CoV-2 (virus that causes COVID-19)
Respiratory syncytial virus
Human Immunodeficiency virus type 1 (HIV-1)*
Duck Hepatitis B Virus (Surrogate for Human Hepatitis B Virus)*
Bovine Viral Diarrhea Virus (Surrogate for Human Hepatitis C Virus)*
Avian Influenza A (H3N2) Reassortant virus
Influenza A (H1N1) virus
Herpes simplex virus type 1
Herpes simplex virus type 2

DISINFECTANT: 2-minute contact time
DÉSINFECTANT: Temps de contact de 2 minutes
Vancomycin Resistant *Enterococcus faecalis* (VRE)
Methicillin Resistant *Staphylococcus aureus* (MRSA)
Salmonella enterica
Trichophyton interdigitale

DISINFECTANT: 3-minute contact time
DÉSINFECTANT: Temps de contact de 3 minutes
Mycobacterium bovis BCG
Pseudomonas aeruginosa
Staphylococcus aureus
Poliovirus Type 1

DIRECTIONS FOR USE

This product is not to be used as a terminal sterilant/ high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

TO START THE FEED THROUGH THE LID:

1. Unseal entire lid from the canister and remove the protective seal from the canister.
 2. Pull up a corner of the first wipe from the center of the roll and push 1 inch through the hole in the lid. DO NOT PUSH FINGER THROUGH THE HOLE.
 3. Replace lid tightly on the canister.
 4. Lift the flip-top and pull out the wipes one at a time.
 5. Ensure lid remains closed between uses to prevent moisture loss.
- This container may be refilled by the consumer with tb Minuteman NEX GEN Disinfectant Wipes Refill Pouch. Refer to further instructions on the label of the refill pouch.

TO SANITIZE HARD, NONPOROUS, NONFOOD CONTACT SURFACES**

Pre-clean surfaces using a wipe prior to sanitization. Discard. Using tb Minuteman NEX GEN Disinfectant Wipes, thoroughly wet the surface and allow to remain wet for 30 seconds. Then wipe the surface dry.

TO DISINFECT HARD, NONPOROUS SURFACES**

Pre-clean heavily soiled surfaces using a wipe. Discard. Using tb Minuteman NEX GEN Disinfectant Wipes, thoroughly wet the surface and allow to remain wet for the appropriate time indicated for the purpose intended. Use additional wipes if needed to ensure the correct wet contact time. For broad spectrum disinfection, allow the surface to remain wet for 3 minutes. Then wipe the surface dry. Do not contaminate food during use and storage.

****A potable water rinse is recommended after sanitizing or disinfecting if:**
- any surfaces may come into direct contact with children at the mouthing stage, or
- any surfaces may come into contact with food or drinking water.

*This product effectively kills Human Immunodeficiency virus 1 (HIV-1), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV) in settings where the viruses would be expected to be encountered such as settings during contamination by blood or body fluid is likely.

Special Decontamination Procedures:

Personal protection: Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e. barrier protection), such as disposable gloves, gowns, masks or eye coverings;

Cleaning: Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product; and

Infectious materials disposal: Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

#Cleaning Care for Dental Chairs
Call 1-800-268-8633 for more information
www.maxill.com/ca/chaircare

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