For Caries Arrest

Description:

5 mL amber dropper bottle of SilverSense SDF solution. Dosage of 0.03 mL per drop. Product forms insoluble precipitates with calcium or phosphate in the dentinal tubules to block nerve impulses. Arrests the continuation of caries.

Uses:

• Helps with the treatment of dentinal sensitivity • Anti-caries • Prevents, fights and/or protects against cavities or caries • Arrests the progress of an already formed cavity in primary teeth • Arrests the continuation of a cavity that has already formed in primary teeth • Helps arrest the progress of an already formed cavity in permanent teeth • Helps arrest the continuation of a cavity that has already formed in permanent teeth. Helps to temporarily reduce (painful) tooth sensitivity due to air exposure in adults.

Contraindications:

1. Patients with ulcerative gingivitis or stomatitis or known sensitivity to silver or other heavy-metal ions.
2. Patients having had full mouth gingivectomies.
3. Patients showing abnormal skin sensitization in daily circumstances.

Directions For Use:

Adults and Children 3 years and older: Use 1-2 drops per patient (0.03 ml per drop). For professional use only.

1. Isolate the affected tooth.
2. Clean and dry the affected tooth surface.
3. Dispense 1-2 drops of solution into a disposable dappen dish (treats up to 5 sites per patient). Transfer material directly to the tooth surface with the supplied Benda® Micro applicator. Use care to avoid getting the solution onto soft tissue or spilling on counters or clothing.
4. Allow to air dry for 60 seconds, do not rinse. If needed, one or two re-applications may be administered at intervals of one week.
5. Do not eat or drink 30-60 minutes after treatment.

Note:

When possible, we suggest removing soft infected dentin with an excavator, or similar instrument, to reduce the potential for staining.

Warnings:

SilverSense SDF is intended for local application only. Not for ingestion. Protect the patient's eyes. Use caution to avoid contact with skin or clothing. In the event of exposure to eyes or skin, flush the area copiously with water and immediately seek medical consultation. Treatment can cause lesions (cavities) in teeth to be permanently stained dark brown or black. Do not use in patients with oral soft tissue ulcerations.

Adverse Reactions:

1. SilverSense SDF does not normally stain enamel or burnished dentin. Advise patients that soft dentin or margins of composite restorations may be stained. Staining may be reversed by gentle polishing with a tincture of iodine (weak iodine solution).
2. Advise patients that air-drying and product application can cause momentary transient pain to hypersensitive areas.
3. Minimize product contact with gingiva and mucous membrane by using recommended amounts and careful application, as transient but rare irritation of the gingiva has been reported. Apply petroleum jelly or cocoa butter, or use cotton rolls, to protect the gingival tissues. Alternatively, a rubber dam can be used to isolate the area.

Precautions For Handling:

1. Store in original packaging in a cool, dark place. Do not freeze or expose to extreme heat.
2. Replace cap immediately after use. Use as soon as dispensed.
3. SilverSense SDF will stain skin, clothes, countertops, floors and instruments brown or black. Refer to the following for stain removal:
 - Skin; wash immediately with water, soap, ammonia or iodine tincture and then rinse thoroughly with water. Do not use excessive methods in an attempt to remove difficult stains from skin as the stains will eventually fade.
 - Clothing/Countertops/Floors/Instruments; use the same procedure as with stained skin. Difficult stains may be treated with sodium hypochlorite.
 - Use caution to avoid transferring the material from gloved hands to other surfaces.
4. After SilverSense SDF is dispensed, dispose of or be sure to wash or thoroughly wipe the container clean immediately after use
5. SilverSense SDF has a blue tint to increase visibility in dappen dish.
6. When dosing, ensure SilverSense SDF is kept separate and not confused with other colored liquids, such as etchants and tooth conditioners.

Medicinal Ingredients:

Silver 25.5% w/v (Silver Diamine Fluoride), Fluoride 4.5% w/v (Silver Diamine Fluoride)

Non-Medicinal Ingredients:

Distilled Water, Fast Green FCF, Ammonium hydroxide.

Visit centrixdental.com/IFU-SDS For Safety Data Sheet

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