maxill®

Tartar & Stain Remover

Cleans away food & tobacco stains and calculus from dentures, bridges & orthodontic appliances and removes rust, calcium hydroxide and hard water scale.

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Tartar and Stain Remover

WARNING

Harmful if swallowed.
Causes serious eye irritation.

PREVENTION
Wash hands after handling. Do not eat, drink or smoke when using this product. Wear eye/face protection.

RESPONSE
IF SWALLOWED: Get medical help. Rinse mouth.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical help.

STORAGE
Keep out of reach of children.

DISPOSAL
Dispose of contents, wastes and residues in accordance with local/regional/national regulations.

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Voir la FDS

Directions: Tartar & Stain Remover (TSR) is ready-to-use. Pour an adequate amount of TSR into a glass beaker suitable to a level that will totally immerse the items to be cleaned. Items may be cleaned by soaking and brushing. If ultrasonic cleaning is desired, cover the beaker and place it into an ultrasonic unit with 3 cm of the beaker sitting in the carrier bath. Ensure carrier bath will not enter the beaker. Sonic clean for 5 to 10 minutes. Immediately thoroughly rinse items with water. Mouth items, such as dentures, etc., should be neutralized in a solution of 1 teaspoon of sodium bicarbonate diluted in a cup of water. Rinse item again. Neutralize waste or unwanted solution using ample amounts of sodium bicarbonate to bring pH above 7.

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+D89260705600100C

maxill inc.
www.maxill.com 1-800-268-8633

Manufactured for:
USA: maxill inc.
Cortland, OH
44410 USA

Canada: maxill inc.
St. Thomas, ON
N5R 6C8
Canada

Made in Canada with domestic and imported parts.

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