DOVE®

maxill

Disposable Valves Tailpiece Options



HVE and SE Valve Tailpieces are designed to securely connect to the HVE or SE tubing to accommodate DOVE Disposable Valves. These tailpieces are simple to install and are compatible with dental chair hangers from all of the most popular manufacturers.

Step 1



Locate the HVE or SE tailpiece to be removed at the end of the vacuum tubing.

Step 2



Cut off the existing tailpiece, approximately two inches below its current position.

Step 3



Push the new DOVE® tailpiece on the existing vacuum tubing, using dish soap or natural lubricant as needed.

Step 4



Install and use the disposable DOVE® HVE or SE valve.



Standard HVE Tailpiece

Item 53085HS

Recommended for the following Dental Chair Hangers:

a-dec · Belmont · forest dental · MIDMARK · Ritter THE DENTAL EXPERTS · TPC



Wide HVE Tailpiece

Item 53085HW

Recommended for the following Dental Chair Hangers:

DentalEZ INTEGRATED SOLUTIONS · Royal · Proma · marus · Pelton & Crane · ENGLE DENTAL SYSTEMS



SE Valve Tailpiece

Item 53084S

DOVE SE Valves fit most manufacturers tailpiece adapters sold in North America. We offer the SE Valve Tailpiece to accommodate unique fits.