OHIO FORGE®

POLISHERS

Directions for Use

COMPOSITE 1-STEP UNIVERSAL POLISHING SYSTEM

1-step diamond polishing system for smoothing and polishing dental composites to the ultimate shine.

Indications:
Ultimate single-step polishing of all composite restorations.

Directions for use:
Trim and smooth restoration with **diamond** or **carbide finishing bur** to maximum smoothness.
1) **Polish** beginning with firm pressure and higher speed. Maximum rpm: 15,000. Wipe restoration with cotton. Continue polishing with reduced pressure and speed with light touch to create desired shine.
2) For best results, polish wet.
Composition: All polishers contain a synthetic rubber matrix infused with diamond particles in various sizes and pigments (mainly Titanium Dioxide). The shanks are made of stainless steel and are surface refined with gold flashing.

COMPOSITE 2-STEP POLISHING SYSTEM

2-step diamond polishing system for smoothing and polishing dental composites to a high-luster.

Indications:
High-luster polishing of all composite restorations.

Directions for use:
1) Reduce and smooth with **medium grit (pink)** polisher.
2) Finish to high-luster with **fine grit (white)** polisher.
IMPORTANT: Clean restoration surface with gauze between polishing steps to remove polish residues.
Polish wet at 8,000 rpm. Maximum rpm: 15,000.
Composition: All polishers contain a synthetic rubber matrix infusedwith diamond particles in various sizes and pigments (mainly Titanium Dioxide). The shanks are made of stainless steel.

COMPOSITE 1-STEP SINGLE-PATIENT POLISHING SYSTEM

1-step diamond polishing system for smoothing and polishing dental composites to a high-gloss.

Indications:
High-gloss single-step polishing of all composite restorations.

Directions for use:
1) **Polish** inished restoration until desired shine is achieved.
Polish wet at 8,000 rpm. Maximum rpm: 15,000.
Composition: All C1S polishers contain a synthetic rubber matrix infused with diamond particles in various sizes and pigments (mainly Titanium Dioxide). The shanks are made of high performance polymer.

OHIO FORGE®

POLISHERS

PORCELAIN & METAL 3-STEP POLISHING SYSTEM

3-step polishing system for smoothing and polishing dental porcelain and metals to a high-luster.

Indications:
High-luster polishing of all porcelain & metal restorations.

Directions for use:
1) Reduce and trim with **coarse grit (green)** polisher.
2) Smooth with **medium grit (pink)** polisher.
3) Finish to high-luster with **fine grit (white)** polisher.
Polish wet or dry at 8,000 rpm. Maximum rpm: 15,000.
Composition: All polishers contain synthetic rubber, diamond grain in various sizes and pigments (mainly Titanium Dioxide). The shanks are made of stainless steel.

PORCELAIN & METAL 2-STEP SINGLE-PATIENT POLISHING SYSTEM

2-step polishing system for smoothing and polishing dental porcelain and metals to a high-gloss.

Indications:
High-gloss polishing of all porcelain & metal restorations.

Directions for use:
1) Reduce and smooth with **medium grit (pink)** polisher.
2) Finish to final gloss with **fine grit (yellow)** polisher.
Polish wet or dry at 8,000 rpm. Maximum rpm: 15,000.
Composition: All polishers contain a silicone matrix infused with silicon carbide particles in various sizes and pigments (mainly Titanium Dioxide). The shanks are made of high quality polymer.

Warnings:
The product should not be used if a patient is known to be allergic to any of the ingredients of the polishing instruments.
Avoid inhaling grinding dust.
Do not exceed maximum rotary speed.

Storage:
Storage temperature 2-28°C/36-82°F.
Do not expose to direct sunlight.
Keep out of the reach of children!
For use in dentistry only!

®maxill inc.
www.maxill.com 1-800-268-8633

Manufactured for:
USA: maxill inc.
Cortland, OH
44410 USA

Canada: maxill inc.
St. Thomas, ON
N5R 6C8
Canada

REV.2023/03