




maxill u-test PCD

Process Challenge Device for Sterilization Monitoring

A PCD is a key element in the quality assurance testing of dental office sterilizers. It is used to monitor the performance of the sterilization process. The PCD simulates an equal or greater challenge than the most difficult instrument/item routinely processed in a sterilization cycle.

Using the **maxill u-test PCD** adds a more rigorous approach to testing and follows best practice within the risk management of the theory of routine practice principles of infection prevention and control.



Blue



Purple



Pink



Teal



Simulates the most difficult instrument within the sterilization load



Compact design saves space in the sterilizer



Manufactured with lightweight, durable aluminum and lasts for an unlimited number of cycles without a change in performance.



Easy-to-use, consistent, reproducible challenge

PRODUCT AVAILABILITY

Blue PCD	#81155
Purple PCD	#81156
Pink PCD	#81157
Teal PCD	#81158

ADDITIONAL PRODUCT INFO



Learn more by watching product information videos online at www.maxill.com

Made in China

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of laboratory testing.