



## maxill xacto stripz sawz

### Interproximal Reduction Strips

The xacto stripz sawz interproximal reduction (IPR) system is indicated mainly for orthodontic stripping (technique for creating interproximal space by wearing down the enamel via abrasion) and can also be used for interproximal finishing and for removal of excess cement and restorative materials. It has a sequence of five strips - one saw and four abrasive strips ranging from extra fine to coarse.

The xacto stripz sawz strips can be used on their own, or inserted into a reciprocating handpiece or optional xacto stripz sawz handle.



### Product Availability

#### Serrated

0.05 mm - SAW ..... #81450  
(10 strips per pack)

#### Right-Sided

0.05 mm - EXTRA FINE ..... #81451  
0.10 mm - FINE ..... #81452  
0.15 mm - MEDIUM ..... #81453  
0.20 mm - COARSE ..... #81454  
(10 strips per pack)

#### Left-Sided

0.05 mm - EXTRA FINE ..... #81455  
0.10 mm - FINE ..... #81456  
0.15 mm - MEDIUM ..... #81457  
0.20 mm - COARSE ..... #81458  
(10 strips per pack)

#### Double-Sided

0.10 mm - EXTRA FINE ..... #81460  
0.20 mm - FINE ..... #81461  
0.30 mm - MEDIUM ..... #81462  
0.40 mm - COARSE ..... #81463  
Assorted Kit (2 of each grit + 2 saws) ..... #81466  
(10 strips per pack)

#### Handle

xacto stripz sawz Handle  
#81449



Item #	#81450	#81451(Right) or #81455 (Left)	#81460	#81452(Right) or #81456 (Left)	#81461	#81453(Right) or #81457 (Left)	#81462	#81454(Right) or #81458 (Left)	#81463
Interproximal Space Right or Left Sided	0.05 mm	0.05 mm	—	0.10 mm	—	0.15 mm	—	0.20 mm	—
Interproximal Space Double Sided		—	0.10 mm	—	0.20 mm	—	0.30 mm	—	0.40 mm
Width	4 mm	4 mm		4 mm		4 mm		4 mm	
Length	15 mm	15 mm		15 mm		15 mm		15 mm	
Grit Size	Saw (no abrasive)	Extra Fine		Fine		Medium		Coarse	

Made in Brazil

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied.

MP-REF-81450B  
Rev.0: 2025-05-23