


maxill Plus

EARLOOP STYLE PROCEDURAL MASKS
WITH ANTI-FOG VISOR

maxill Plus Earloop Style Procedural Masks with Anti-fog Visor feature a 4-ply construction and provide Level 3* protection. A soft inner layer, anti-fog foam strip, and flexible nosepiece ensure maximum comfort and a secure fit. These masks are not made with natural rubber latex.



Availability

- Blue#53041
 - Pink#53042
 - Green#53043
 - Lavender#53044
 - White#53045
- (25 masks per box)

***ASTM LEVEL 3** Complies with ASTM F2100 Requirements:

- Fluid Resistance:** 160 mm Hg
- Bacterial Filtration Efficiency:** BFE ≥ 98%
- Particle Filtration Efficiency:** PFE ≥ 98%
- Breathability - Delta P:** < 6.0 mm H₂O/cm²
- Flame Spread:** Class 1

How to Wear an Earloop Mask



1 Use your fingers to open the ear-loops. Position the nose clip at the top. Place the mask over your nose and mouth.



2 Position the ear-loops around your ears. Pull the mask from the top and bottom to unfold the pleats of the mask.



3 Form the nose clip over the bridge of your nose to make it fit snugly.

Made in China

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-53024B
Rev.1: 2025-10-09