

SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - Product/ Company Identification

Product Trade Name:
Tempit® Ultra-F - Light-Cure Temporary Restorative Paste with Fluoride
Tempit® Ultra - Light-Cure Temporary Restorative Paste
Synonyms: Light-Cure Temporary Restorative with/without Fluoride Paste
Recommended Use: Dental Material
Part (Item) Number: 310064, 310065
Manufacturer: CENTRIX INC, 770 RIVER ROAD, SHELTON, CT 06484 U.S.A.
Emergency Telephone Number: (800) 235 - 5862

SECTION II – HAZARDS IDENTIFICATION

For potential health hazards: SEE Section XI. 
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents:
Tempit® Ultra-F - Light-Cure Temporary Restorative Paste with Fluoride
Tempit® Ultra - Light-Cure Temporary Restorative Paste

Hazardous Ingredients:

Ingredient	HAZARDOUS DATA	Hazardous Data Legend
Sodium Fluoride Solid	Health: 3	4 - Extreme
CAS: 7681-49-4	Flammability: 0	3 - High
UN Designation: UN1690	Reactivity: 1	2 - Moderate
Class: 6.1. Toxic		1 - Slight
Packing Group: III		0 - Insignificant
Exact % content withheld due to trade secret restrictions.		

SECTION IV - FIRST AID MEASURES

SKIN: Wash with soap and water.
INGESTION: Induce vomiting if large amount of material is swallowed.
INHALATION: No harm anticipated.
EYES: Flush with copious amounts of water.

SECTION V - FIRE-FIGHTING MEASURES

SUITABLE EXTINGUISHING MEDIA: Water contact acceptable.
FORBIDDEN EXTINGUISHING MEDIA: None.
UNUSUAL FIRE AND EXPLOSION HAZARDS: None specific.
SPECIAL PROTECTION EQUIPMENT: None.

SECTION VI - ACCIDENTAL RELEASE MEASURES

PERSONAL PRECAUTIONS: Avoid contact of skin or eyes. Wash hands with soap and water if contact occurs.
ENVIRONMENTAL PRECAUTIONS: Dispose of in accordance to all federal, state, local regulations.
RECLAIMING METHODS: Wipe up with soap, water and paper towels.

SECTION VII - HANDLING AND STORAGE

HANDLING: None. STORAGE: Store at temperatures not exceeding 75°F (24°C).
CONTAINERS: Avoid exposure to moisture and light.
INCOMPATIBILITIES: None specific. OTHER PRECAUTIONS: None.

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

GENERAL MEASURES: None.
SKIN PROTECTION: Avoid prolonged contact skin; wash with soap and water.
EYE PROTECTION: Goggles/Eyeshields.
HANDS PROTECTION: Avoid prolonged contact skin; wash with soap and water.
LIMITS: Respiratory protection: Not required
OTHER EQUIPMENT:

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE: Tooth shade paste
ODOR: Negligible odor
BOILING POINT: N/A
EXPLOSION LIMITS: N/A
WATER SOLUBILITY: Virtually insoluble
EVAPORATION RATE: Not established
FLASH POINT: N/A

SECTION X - STABILITY AND REACTIVITY.

STABILITY: Generally stable.
CONDITIONS TO AVOID: Light exposure, heat.
MATERIALS TO AVOID: N/A
HAZARDOUS DECOMPOSITION PRODUCTS: N/A - unless exposed to high temperatures and/or light and only when large quantities of the material are involved.

SECTION XI - TOXICOLOGICAL INFORMATION.

POTENTIAL ACUTE HEALTH EFFECTS: Possibility of skin irritation or sensitization if not washed off promptly.
POTENTIAL CHRONIC HEALTH EFFECTS: Possibility of skin sensitization or irritation upon prolonged or repeated exposure.
SIGNS AND SYMPTOMS OF EXPOSURE: Skin/eye irritation.
MEDICAL CONDITIONS GENERALLY AGGRAVATE: None known.
CARCINOGENIC OR POTENTIAL CARCINOGENIC: No ingredient listed as a carcinogen.
OTHER POTENTIAL HEALTH EFFECTS: None known.

SECTION XII - ECOLOGICAL INFORMATION.

No ecological damage or impact anticipated when handled as indicated.

SECTION XIII - DISPOSAL CONSIDERATIONS.

WASTE DISPOSAL: Do not dispose into sink. Dispose in accordance with all federal, state, and local regulations.

SECTION XIV - TRANSPORT INFORMATION.

No hazardous goods within the meaning of transport regulations.

SECTION XV - REGULATORY INFORMATION

U.S. FDA DENTAL DEVICE: Class II product
EC MEDICAL DEVICE DIRECTIVE: Class IIa product

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION: The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility

centrix®
Making Dentistry Easier.™



CENTRIX INCORPORATED

770 River Road, Shelton, CT 06484 USA
(203) 929-5582 • (800) 235-5862 • FAX (203) 929-6804
www.centrixdental.com



Dental-Contact Vertriebs OHG, Im Spiet 101, DE 26506 Norden

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia