

maxill

# tb Minuteman NEX GEN DISINFECTANT WIPES



## MYCOBACTERICIDAL • VIRUCIDAL<sup>1</sup> • BACTERICIDAL • FUNGICIDAL

tb Minuteman NEX GEN DISINFECTANT WIPES have a unique formulation of four active ingredients that work synergistically to create an unparalleled broad spectrum disinfectant ideal for daily use in dental offices. tb Minuteman NEX GEN DISINFECTANT WIPES can be used on a multitude of hard surfaces to inactivate a broad spectrum of pathogenic microorganisms.

- Broad Spectrum Bactericide, Virucide, Tuberculocide and Fungicide
- Meets CDC guidelines for intermediate level surface disinfection
- Cleaner and disinfectant in one product; save money, time and storage space
- Compatible with a wide range of hard, non-porous surfaces
- Low lint, 6" x 7" heavy duty, non-woven, pre-saturated wipes
- Packaged in recyclable HDPE tubs with self-closing dispenser lids
- Convenient refill rolls are available to reduce waste and save you money




### Heavy Duty Cloth

Heavy-duty, non-woven cloth design provides uniform liquid release on surfaces and prevents bunching or tearing in the dispenser.



### One-Stop Shop

tb Minuteman NEX GEN DISINFECTANT is a versatile disinfectant solution that can be used on a wide variety of surfaces throughout the establishment.



### Mountable Rack

Available wall mountable tub rack (sold separately) facilitates increased compliance by making it easier to have a disinfectant wipe within reach.



### Broad Spectrum Disinfection

The wipes are not only effective at cleaning, they also provide broad spectrum disinfection in efficient contact times. Soiled surfaces require a preclean step.



### Reduce Waste

Convenient refill rolls sold in plastic pouches reduce waste and save money. Our pouches are triple-laminated to help maintain stability until you are ready to use the roll.



### Calendar Stickers

Easy-to-read calendar is provided with every refill pouch to easily keep track of when product was opened, making inventory management easier.

Made in Canada with domestic and imported parts.

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-61160B  
Rev.3: 2023-03-06

maxill inc., St. Thomas, ON N5R 6C8  
[www.maxill.com](http://www.maxill.com) 1-800-268-8633

## † Contact Times

### **SANITIZER:**

#### **30 seconds contact time**

*Staphylococcus aureus*  
*Klebsiella pneumoniae*

### **DISINFECTANT:**

#### **1-minute contact time**

Rotavirus  
Human Coronavirus  
SARS-CoV-2 (virus that causes COVID-19)<sup>†</sup>  
Human Immunodeficiency virus type 1 (HIV-1)\*  
Duck Hepatitis B Virus (Surrogate for Human Hepatitis B Virus)\*  
Bovine Viral Diarrhea Virus (Surrogate for Human Hepatitis C Virus)\*  
Avian Influenza A (H3N2) Reassortant virus  
Herpes simplex virus type 1  
Herpes simplex virus type 2

#### **2-minute contact time**

Vancomycin Resistant *Enterococcus faecalis* (VRE)  
Methicillin Resistant *Staphylococcus aureus* (MRSA)  
*Salmonella enterica*  
*Trichophyton interdigitale*

#### **3-minute contact time**

*Mycobacterium bovis* BCG  
*Pseudomonas aeruginosa*  
*Staphylococcus aureus*  
Poliovirus Type 1



Designed with **You** in mind!

### **Product Availability**

|                       |        |
|-----------------------|--------|
| 160 wipes tub         | #61160 |
| 160 wipes refill roll | #61190 |

Visit us online for more information,  
including SDSs and pricing:

[www.maxill.com/nexgen](http://www.maxill.com/nexgen)

\*This product effectively kills Human Immunodeficiency virus 1 (HIV-1), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV) on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HIV-1, HBV, and HCV.

Made in Canada with domestic and imported parts.

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-61160B  
Rev.3: 2023-03-06

maxill inc., St. Thomas, ON N5R 6C8  
[www.maxill.com](http://www.maxill.com) 1-800-268-8633