

Ultrasonic Bath:

Cleaning, Solution Changing & Test Log

Dental Office:	Ultrasonic Bath#:
IPAC Rep:	(each bath requires its own logging form)

Date	Cleaning & Disinfecting of Bath	Solution Filled in am	Solution Emptied in pm	*Extra Solution Change	Strip Test	Strip Test Results Pass or Fail	Initials

^{*}Best practice for optimal efficacy indicates changing the enzymatic solution midday, especially for large practices processing a high volume and/or soiled surgical instruments.