ITEM 60087



For cleaning and control of microbial contaminants in dental unit waterlines

Purifies water for dental delivery units, filtration equipment, independent water bottles, scalers and other dental process-water applications.

KEEP OUT OF REACH OF CHILDREN WARNING

See label for additional precautionary statements

Active Ingredient:

 Sodium Dichloroisocyanurate
 16.88%

 Other Ingredients:
 83.12%

 Total:
 100.00%

EPA Reg. No. 71847-4-90841 EPA Establishment No. 71847-IRL-001

Net Contents: 50, 49 mg Tablets



DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Read the entire label and use strictly in accordance with precautionary statements and directions

This product is effective against microbial contaminants

TREATING DENTAL PROCESS-WATER AT ROOM TEMPERATURE

To Use: Remove one maxill H₂O DUWL TABS tablet from packaging and quickly insert into 1/5 gallon or up to 2 quarts of water. Ensure tablet is fully dissolved before use.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage and disposal.

Pesticide Storage: Store in a cool, dry area away from sunlight and heat and in an area inaccessible to children.

Pesticide Disposal: Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer.

Container Handling: Nonrefillable packet. Do not reuse or refill this packet. Offer for recycling, if available.

Emergency Handling: In case of contamination or decomposition; do not reseal container. If possible isolate container in open and well-ventilated area. Flood with large volumes of water. Dispose of contaminated material in an approved landfill area.

FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

If swallowed: Call a poison control center, or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person.

If inhaled: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. Call a poison control center or doctor for further treatment advice.

If on skin or clothing: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15–20 minutes. Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

YOU MAY ALSO CONTACT 1-800-222-1222 FOR EMERGENCY MEDICAL TREATMENT INFORMATION.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric layage

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

WARNING: Causes substantial but temporary eye injury. Harmful if swallowed.

Do not get in eyes or on clothing. Wear protective eyewear (goggles with side-shields, safety glasses, or face shield). Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS

This pesticide is toxic to fish and aquatic organisms.

PHYSICAL OR CHEMICAL HAZARDS

Contains strong oxidizing agent. Use only clean dry utensils. Mix only into water. Contamination with moisture, dirt, organic matter or other chemicals or any other foreign matter may start a chemical reaction with generation of heat, liberation of hazardous gases and possible generation of fire and explosion. Avoid any contact with flaming or burning materials such as a lighted cigarette. Do not use this product in any chlorinating device which has been used with any inorganic or unstabilized chlorinating compounds (e.g., calcium hypochlorite). Such use may cause fire or explosion.

Manufactured for:

*maxill inc.
500 W. Main St., Corlland, OH, 44410
1-855-462-9455 www.maxill.com

RFV 2019/10

Made in Ireland