MAXII Procedural

maxill inc., St. Thomas, ON N5R 6C8 www.maxill.com 1-800-268-8633

maxill chrömamax

High Precision Chromatic Alginate

maxill chrömamax High Precision Chromatic Alginate features chromatic indication of the working phases for easy and precise results.

Working Time (23° C) Including mixing time	1m 25s
Time in the mouth	30s
Dimensional Stability	7 Days
Flavour	Spearmint
Final Colour	Blue

- Long Dimensional Stability
- Easy-to-use



- High Elasticity & Resistant to Tear
- Long, 5 year Expiry Date*



Fast Setting Time

Chromatic Phase Indicator



Product Availability

chrömamax 450 g bag	#53775
chrömamax Storage Canister	#53775T

Made in Italy

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.