

## **Zymax Enzymatic Cleaning Solution** is an eco-friendly multipurpose cleaner.

A blend of enzymes in a low foaming solution, Zymax Enzymatic Cleaning solution loosens and removes dried debris, blood, protein and mucin deposits. Zymax is also used as a pre-cleaner in the maxill H<sub>2</sub>O Dental Unit Waterline Treatment and Monitoring System.

- Non Corrosive
- Non Caustic
- ► No Phosphates
- Low Foaming

## Use Zymax to safely clean:\*

- Plastics
- Rubber
- Stainless steel
- Fabrics soiled with organic debris

## **Product Availability**

**Zymax** - 32 fl. oz. Spray Bottle (Ready-to-Use) . . . . #60014 **Zymax 128** - 1.06 gt. Bottle (128x Concentrate) . . . #60035

128 L

Zymax 128 concentrate will make up to 128 L of working solution







946 mL (32 FL OZ)

\*Zymax will safely clean synthetic and natural fabrics that are soiled with organic debris. Prolonged contact with a Zymax working solution with certain materials may cause discoloring or corrosion. Prior to using a Zymax working solution, the user should first perform compatibility testing. Monitor contact times according to the label.

Made in Canada with domestic and imported parts

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-60013BU Rev.3: 2021-03-05