





Zymax QS Enzymatic Cleaning Foam is a multi-enzyme RTU detergent. The proprietary blend of enzymes is formulated to be especially effective in the removal of organic material found on instruments.

Zymax QS is ideal for saturating instruments when soaking is not an immediate option. The spray foam coats the instrument and reaches hard-to-reach places to help ensure organic matter does not harden and complicate the cleaning/reprocessing process. Zymax QS Enzymatic Cleaning Foam is a complementary addition to instrument re-processing as a pre-cleaning step.



Use Zymax QS to safely clean:

- Plastics
- Rubber
- Stainless steel
- Fabrics soiled with organic debris
- Non Corrosive
- No Phosphates
- Non Caustic
- Foaming

Product Availability

Zymax QS - 32 fl oz Spray Bottle (Ready-to-Use) #60025









Made in Canada with domestic and imported parts

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing. Therefore, each health care facility should conduct its own efficacy testing in order to verify the efficacy of the product under its own normal operating conditions to meet the specific requirements of the workplace

MP-REF-60025BU Rev0: 2019-08-20