MAXII Infection Control

maxill inc., Cortland, Ohio 44410 www.maxill.com 1-855-4<u>62-9455</u>



The **Ster-FX BIN** is a germicide tray designed for soaking instruments in a variety of solutions including enzymes, detergents, phenolics, glutaraldehydes and other sterilizing, disinfecting and pre-soaking solutions.

Pair up the Ster-FX BIN with Zymax 128 or Zymax LF Enzymatic Cleaning Solutions to soak instruments before cleaning and sterilization. The Ster-FX BIN can also be used with Ster-FX Chemosterilant & Sporicide to sterilize reusable semi-critical and critical medical devices by immersion.

- Dials on lid can be used to keep track of when items are to be removed and/or when the solution needs to be changed
- Interior tray features a perforated diamond design allowing for quick drain-off and easy transfer of instruments
- > 2 tier tray provides convenient hands free draining
- Holds 1 gallon of liquid and has easy pour corners for disposal of solution
- Hinged, fitted cover helps prevent spills and contain odors
- Inner tray dimensions: 11-1/2" x 6-3/4" x 3"
- Available 1000 mL Ster-FX E-Z Beak Measuring Beaker makes measuring holding solutions quick and easy

Product Availability

Ster-FX BIN	#81108
Ster-FX E-Z BEAK	#81109





USE WITH





