

maxill L1 MEDICAL/DENTAL

EARLOOP STYLE PROCEDURAL MASKS

maxill L1 Medical/Dental Earloop Procedural Masks are not made with natural rubber latex and offer Level 1 protection*.

*ASTM LEVEL 1

Complies with ASTM F2100-19 Requirements:

Fluid Resistance: 80 mm Hg

Bacterial Filtration Efficiency: BFE $\geq 95\%$ **Particle Filtration Efficiency:** PFE $\geq 95\%$ **Breathability - Delta P:** $< 4.0 \text{ mm H}_2\text{O/cm}^2$

Flame Spread: Class 1

Availability

Blue.....#53017



How to Wear an Earloop Mask



Use your fingers to open the ear-loops. Position the nose clip at the top. Place the mask over your nose and mouth.



Position the ear-loops around your ears. Pull the mask from the top and bottom to unfold the pleats of the mask.



Form the nose clip over the bridge of your nose to make it fit snugly.

Made in China