



maxill Silken^{4P}

4-PLY EARLOOP STYLE PROCEDURAL MASKS

maxill Silken 4-Ply Earloop Procedural Masks are not made with natural rubber latex and offer Level 3 protection*.

*ASTM LEVEL 3

Complies with ASTM F2100-19 Requirements:

Fluid Resistance: 160 mm Hg

Bacterial Filtration Efficiency: BFE ≥ 98%

Particle Filtration Efficiency: PFE ≥ 98%

Breathability - Delta P: < 5.0 mm H₂O/cm²

Flame Spread: Class 1

Availability

Blue.....	53056
Pink.....	53057
Green.....	53058
Lavender.....	53059
White.....	53060

(50 masks per box)



How to Wear an Earloop Mask



1 Use your fingers to open the ear-loops. Position the nose clip at the top. Place the mask over your nose and mouth.



2 Position the ear-loops around your ears. Pull the mask from the top and bottom to unfold the pleats of the mask.



3 Form the nose clip over the bridge of your nose to make it fit snugly.