



## **4-PLY EARLOOP STYLE PROCEDURAL MASKS**

maxill Silken 4-Ply Earloop Procedural Masks are not made with natural rubber latex and offer Level 3 protection\*.

## \*ASTM LEVEL 3

Complies with ASTM F2100-19 Requirements:

Fluid Resistance: 160 mm Hg

**Bacterial Filtration Efficiency:** BFE  $\ge 98\%$ **Particle Filtration Efficiency:** PFE  $\ge 98\%$ **Breathability - Delta P:**  $< 5.0 \text{ mm H}_2\text{O/cm}^2$ 

Flame Spread: Class 1

## **Availability**

Blue	53056
Pink	53057
Green	53058
Lavender	53059
White	53060

(50 masks per box)



## How to Wear an Earloop Mask



Use your fingers to open the ear-loops. Position the nose clip at the top. Place the mask over your nose and mouth.



Position the ear-loops around your ears. Pull the mask from the top and bottom to unfold the pleats of the mask.



Form the nose clip over the bridge of your nose to make it fit snugly.