



CHAIR CLEANER & POLISH

SCHIFF'S CHAIR CLEANER & POLISH may be used on non-porous environmental surfaces and smooth upholstery for a streak free shine. Its formulation is gentle and safe for use on leather & vinyl upholsteries, stainless steel appliances and fixtures, Formica or stone-type solid counters and tables, electronics, plastics and more.

SCHIFF'S CHAIR CLEANER & POLISH does not contain mineral oil. is low in VOCs and its non-flammable formula helps to prevent static build-up.

CLEANS & POLISHES:

- Dental Chairs
- Work Surfaces





PRODUCT AVAILABILITY

946 mL Spray Bottle#60021 **160 Wipes Tub**.....#60031



Cleaning Care for Dental Chairs

In addition to disinfecting, a cleaning protocol is required by the chair manufacturer. To learn more about how SCHIFF'S CHAIR CLEANER & POLISH can be used to help prolong the life of your dental chair upholstery, call 1-855-462-9455 or visit www.maxill.com/us/chaircare

Made in Canada.