


















NITRILE GLOVES	Glove Name	Colour	Packaging (10 boxes per case)	Total Length (mm)	Width (mm)						Thickness mm/mil	
					XS	S	M	L	XL	XXL	Palm	Finger
	maxill 2020		200 per box	230	N/A	85 ± 5	95 ± 5	105 ± 5	115 ± 5	N/A	0.07/2.8	0.08/3.1
	maxill BK5		100 per box	240	75 ± 5	85 ± 5	95 ± 5	105 ± 5	115 ± 5	N/A	0.08/3.2	0.10/3.9
	Cobalt Elite		100 per box	240	N/A	85 ± 5	95 ± 5	105 ± 5	115 ± 5	120 ± 5	0.09/3.5	0.17/6.7
	Cobalt 911		100 per box	300	N/A	85 ± 5	95 ± 5	105 ± 5	115 ± 5	N/A	0.08/3.2	0.13/5.1
	Coral 200		200 per box	220	70 ± 10	80 ± 10	95 ± 10	N/A	N/A	N/A	0.05/2.0	0.05/2.0
	maxill Nitrile		100 per box	240	N/A	85 ± 5	95 ± 5	105 ± 5	115 ± 5	125 ± 5	0.08/3.1	0.10/3.9
	Nitriclear Aloe AF		250 per box (XL 200)	240	75 ± 5	85 ± 5	95 ± 5	106 ± 5	115 ± 5	N/A	0.05/2.0	0.08/3.1
	nes Earth AF		200 per box	240	75 ± 5	85 ± 5	95 ± 5	106 ± 5	115 ± 5	N/A	0.05/2.0	0.08/3.1
	nes Europa		300 per box (XL 250)	230	70 ± 10	80 ± 10	95 ± 10	110 ± 10	120 ± 10	N/A	0.05/2.0	0.05/2.0
	nes Ceres		200 per box	240	75 ± 5	85 ± 5	95 ± 5	106 ± 5	115 ± 5	N/A	0.06/2.4	0.08/3.1
	nes Pluto		300 per box (XL 250)	230	75 ± 5	85 ± 5	95 ± 5	105 ± 5	115 ± 5	N/A	0.05/2.0	0.06/2.5
	Packard Healthcare		200 per box	240	N/A	85 ± 5	95 ± 5	105 ± 5	N/A	N/A	0.07/2.8	0.10/3.9
	Glove Name	Colour	Packaging	Characteristics								
	Dental Instrument Processing Gloves		100 per box	Puncture Resistance Level 1 NF EN388:2004 §6.3			Dexterity 5 EN420:2003+A:2009 §6.2			Tear Resistance 0 NF EN388:2004 §6.4		
All maxill nitrile gloves are powder free												

NEOPRENE	Glove Name	Colour	Packaging (10 boxes per case)	Total Length (mm)	Width (mm)					Thickness mm/mil		
					XS	S	M	L	XL	Palm	Finger	
	Neoprene EKB Lime		200 per box	240	76 ± 4	86 ± 4	95 ± 4	106 ± 4	115 ± 4	0.07/2.8	0.09/3.5	
All maxill neoprene gloves are powder free												

LATEX	Glove Name	Colour	Packaging (10 boxes per case)	Total Length (mm)	Width (mm)					Thickness mm/mil		
					XS	S	M	L	XL	Palm	Finger	
	maxill PFL1		100 per box	230	70 ± 10	80 ± 10	95 ± 10	111 ± 10	N/A	0.08/3.1	0.08/3.1	
All maxill latex gloves are powder free												

VINYL	Glove Name	Colour	Packaging (10 boxes per case)	Total Length (mm)	Width (mm)					Thickness mm/mil		
					XS	S	M	L	XL	Palm	Finger	
	maxill Vinyl	Colourless	100 per box	240	N/A	85 ± 5	95 ± 5	105 ± 5	115 ± 5	0.11/4.3	0.09/3.5	
	Sensor-maxx		100 per box	230	75 ± 5	85 ± 5	95 ± 5	105 ± 5	N/A	0.08/3.1	0.06/2.4	
	Ultra-maxx		100 per box	230	75 ± 5	85 ± 5	95 ± 5	105 ± 5	115 ± 5	0.08/3.1	0.06/2.4	
All maxill vinyl gloves are powder free												

*While reasonable efforts are made to ensure an accurate representation, glove colours may not be exactly as shown.



Protection Against Fentanyl

Permeation Testing (ASTM D6978-05) for Fentanyl Citrate Injection, demonstrates no breakthrough for up to 240 minutes.



Chemotherapy Tested

Permeation Testing (ASTM D6319) for assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

