


Pair up the Ster-FX BIN with Zymax 128 or Zymax LF Enzymatic Cleaning Solutions to soak instruments before cleaning and sterilization. The Ster-FX BIN can also be used with Ster-FX Chemosterilant & Sporicide to sterilize reusable semi-critical and critical medical devices by immersion.



- Interior tray features a perforated diamond design allowing for quick drain-off and easy transfer of instruments
- 2 tier tray provides convenient hands free draining
- Holds 4 litres of liquid and has easy pour corners for disposal of solution
- Hinged, fitted cover helps prevent spills and contain odours
- ▶ Inner tray dimensions: 11-1/2" x 6-3/4" x 3"
- Available 1000 mL Ster-FX E-Z Beak Measuring Beaker makes measuring holding solutions quick and easy

Product Availability

Ster-FX BIN #81108 Ster-FX E-Z BEAK#81109











