



maxill u•test JUMBO PCD



Process Challenge Device for Sterilization Monitoring

A Process Challenge Device (PCD) is a crucial component in the quality assurance testing of dental office sterilizers. It is designed to help monitor the effectiveness of the sterilization process. The PCD simulates a challenge equal to or greater than the most difficult instrument or item routinely processed during a sterilization cycle.

Using the **maxill u•test JUMBO PCD** provides a more rigorous approach to testing and adheres to best practices in risk management, aligned with the principles of infection prevention and control in routine practice.



Simulates the most difficult instrument within the sterilization load.



Manufactured with lightweight, durable aluminum and lasts for an unlimited number of cycles without a change in performance.



Easy-to-use, consistent, reproducible challenge.



Jumbo size fits larger BIs and CIs that may not fit in the standard size u•test PCD.

Inner diameter: 26 mm

Inner length: 76 mm

PRODUCT AVAILABILITY

PCDs

Blue JUMBO PCD	#81140
Purple JUMBO PCD.....	#81141
Pink JUMBO PCD	#81142
Teal JUMBO PCD.....	#81143

Replacement Caps

Blue JUMBO Replacement Cap.....	#81145
Purple JUMBO Replacement Cap . . .	#81146
Pink JUMBO Replacement Cap.....	#81147
Teal JUMBO Replacement Cap.....	#81148

Made in China

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of laboratory testing.