

# Two Great Options. Same Proven Performance.

## Type 5 Chemical Indicators for Steam Sterilization

maxill offers two excellent Type 5 Chemical Indicator options, ensuring you have reliable choices to meet your sterilization needs.



**u-test SI-LR**



**u-test SI-AB**

Meets ISO 11140-1 Validation Testing	✓	✓
Type 5 Internal Indicator	✓	✓
Integrating Indicator: Time, Temperature and Pressure	✓	✓
Dimensions	72 mm x 14 mm	52 mm x 21 mm

### u-test SI-LR

Type 5 (Class V) Integrating Package Monitoring Indicators

- ✓ Colour Change Bar to Show Results
- ✓ Label-Ready Adhesive Backing
- ✓ Fits Standard or JUMBO PCD

260 pack . . . . . #81125    5,000 pack . . . . . #81121  
 500 pack . . . . . #81117    10,000 pack . . . . . #81127  
 1,000 pack . . . . . #81118  
 2,500 pack . . . . . #81119



### u-test SI-AB

Type 5 (Class V) Migrating Chemical Indicators for Sterilization

- ✓ Window with Migrating Bar to Show Results
- ✓ Fits JUMBO PCD

250 pack . . . . . #81153  
 1,000 pack . . . . . #81154

