


















| NITRILE GLOVES | Glove Name | Color | Packaging (10 boxes per case) | Total Length (mm) | Width (mm) | | | | | | Thickness mm/mil | |
|---|-------------------------------------|---|----------------------------------|---|------------|---------|---------------------------------------|----------|----------|---|------------------|----------|
| | | | | | XS | S | M | L | XL | XXL | Palm | Finger |
| | maxill 2020 |  | 200 per box | 230 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | N/A | 0.07/2.8 | 0.08/3.1 |
| | maxill BK5 |  | 100 per box | 240 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | N/A | 0.08/3.2 | 0.10/3.9 |
| | Cobalt Elite |  | 100 per box | 240 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | 120 ± 5 | 0.09/3.5 | 0.17/6.7 |
| | Cobalt 911 |  | 100 per box | 300 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | N/A | 0.08/3.2 | 0.13/5.1 |
| | Coral 200 |  | 200 per box | 220 | 70 ± 10 | 80 ± 10 | 95 ± 10 | N/A | N/A | N/A | 0.05/2.0 | 0.05/2.0 |
| | maxill Nitrile |  | 100 per box | 240 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | 125 ± 5 | 0.08/3.1 | 0.10/3.9 |
| | Nitriclear Aloe AF |  | 250 per box (XL 200) | 240 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 106 ± 5 | 115 ± 5 | N/A | 0.05/2.0 | 0.08/3.1 |
| | nes Earth AF |  | 200 per box | 240 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 106 ± 5 | 115 ± 5 | N/A | 0.05/2.0 | 0.08/3.1 |
| | nes Europa |  | 300 per box (XL 250) | 230 | 70 ± 10 | 80 ± 10 | 95 ± 10 | 110 ± 10 | 120 ± 10 | N/A | 0.05/2.0 | 0.05/2.0 |
| | nes Ceres |  | 200 per box | 240 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 106 ± 5 | 115 ± 5 | N/A | 0.06/2.4 | 0.08/3.1 |
| | nes Pluto |  | 300 per box (XL 250) | 230 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | N/A | 0.05/2.0 | 0.06/2.5 |
| | Packard Healthcare |  | 200 per box | 240 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | N/A | N/A | 0.07/2.8 | 0.10/3.9 |
| | Glove Name | Color | Packaging | Characteristics | | | | | | | | |
| | Dental Instrument Processing Gloves |  | 100 per box | Puncture Resistance Level 1 NF EN388:2004 §6.3 | | | Dexterity 5 EN420:2003+A:2009 §6.2 | | | Tear Resistance 0 NF EN388:2004 §6.4 | | |
| All maxill nitrile gloves are powder free | | | | | | | | | | | | |

| NEOPRENE | Glove Name | Color | Packaging (10 boxes per case) | Total Length (mm) | Width (mm) | | | | | Thickness mm/mil | | |
|--|-------------------|---|----------------------------------|----------------------|------------|--------|--------|---------|---------|------------------|----------|--|
| | | | | | XS | S | M | L | XL | Palm | Finger | |
| | Neoprene EKB Lime |  | 200 per box | 240 | 76 ± 4 | 86 ± 4 | 95 ± 4 | 106 ± 4 | 115 ± 4 | 0.07/2.8 | 0.09/3.5 | |
| All maxill neoprene gloves are powder free | | | | | | | | | | | | |

| LATEX | Glove Name | Color | Packaging (10 boxes per case) | Total Length (mm) | Width (mm) | | | | | Thickness mm/mil | | |
|---|-------------|---|----------------------------------|----------------------|------------|---------|---------|----------|-----|------------------|----------|--|
| | | | | | XS | S | M | L | XL | Palm | Finger | |
| | maxill PFL1 |  | 100 per box | 230 | 70 ± 10 | 80 ± 10 | 95 ± 10 | 111 ± 10 | N/A | 0.08/3.1 | 0.08/3.1 | |
| All maxill latex gloves are powder free | | | | | | | | | | | | |

| VINYL | Glove Name | Color | Packaging (10 boxes per case) | Total Length (mm) | Width (mm) | | | | | Thickness mm/mil | | |
|---|--------------|---|----------------------------------|----------------------|------------|--------|--------|---------|---------|------------------|----------|--|
| | | | | | XS | S | M | L | XL | Palm | Finger | |
| | maxill Vinyl | Colorless | 100 per box | 240 | N/A | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | 0.11/4.3 | 0.09/3.5 | |
| | Sensor-maxx |  | 100 per box | 230 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 105 ± 5 | N/A | 0.08/3.1 | 0.06/2.4 | |
| | Ultra-maxx |  | 100 per box | 230 | 75 ± 5 | 85 ± 5 | 95 ± 5 | 105 ± 5 | 115 ± 5 | 0.08/3.1 | 0.06/2.4 | |
| All maxill vinyl gloves are powder free | | | | | | | | | | | | |

*While reasonable efforts are made to ensure an accurate representation, glove Colors may not be exactly as shown.



Protection Against Fentanyl

Permeation Testing (ASTM D6978-05) for Fentanyl Citrate Injection, demonstrates no breakthrough for up to 240 minutes.



Chemotherapy Tested

Permeation Testing (ASTM D6978-05) for assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

