




## maxill chrömax

## **High Precision Chromatic Alginate**

maxill chrömamax High Precision Chromatic Alginate features chromatic indication of the working phases for easy and precise results.

Working Time (74° F) Including mixing time	1m 25s
Time in the mouth	30s
Dimensional Stability	7 Days
Flavor	Spearmint
Final Color	Blue



- ► Chromatic Phase Indicator
- ► Long Dimensional Stability
- ► Fast Setting Time
- Easy-to-use

- ► High Elasticity & Resistant to Tear
- ► Long, 5 year Expiry Date\*





## **Product Availability**

**chrömamax 15.9 oz bag**.....#53775 chrömamax Storage Canister ..... #53775T

## Made in Italy

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.